

AutoFuser

ACE Medical Co., Ltd.

3F, Ace Ville B/D 25-1, Bomoon-Dong,
7Ga, Sungbuk-Gu, Seoul 136-087, Korea
Tel: +82-2-921-2967 Fax: +82-2-922-9319

DRUGS ACCEPTABLE FOR AUTOFUSER DISPOSABLE INFUSER

The AutoFuser Disposable Infuser meets the requirements for biocompatibility and drug compatibility. All fluid contact parts are ISO 10993-1 compatible.

The AutoFuser medication reservoir is manufactured using silicone rubber, a material recognized as being highly compatible with most medications. Extensive research exists that establishes silicone rubber as an excellent choice for medical products.

Indwelling intravenous and intra-arterial catheters, percutaneous and subcutaneous infusion lines are just some of the examples of medical devices which contain silicone rubber. Administration of chemotherapy agents, narcotic analgesics, antibiotics, nutrition solutions and blood products can all be conducted safely and reliably with silicone-based medical infusion device such as the AutoFuser.

The following drug list is provided for your convenience to help determine the stability of frequently used medications.

DRUG LIST

PAIN DRUGS	ANTIBIOTICS	CYTOTOXICS	ANTIVIRALS	MISCELLANEOUS DRUGS
Bupivacaine	Amikacin	Bleomycin	Acyclovir	Desferrioxamine
Chloroprocaine	Azlocillin	Carboplatln	Foscamet	Dextrose Anhydrous
Diamorphine	Benzylpenicillin	Cisplatin	Gancicovir	Guanylhydrazone
Droperidol	Carbenicillin	Cyclophosphamide	AZT	Heparin
Fentanyl	Cefamandole	Cytarabine		Metoclopramide
Hydromorphone	Cefonocid	Doxorubicin		NaH CO ₃
Lidocaine	Ceftazidime	Epirubicin		Ondansetron
Meperidine	Ciprofloxacin	5-FU		Salbutamol
Morphine	Erythromycin	Ifosfamide		Sodium Chloride
	Fluconazole	Leucovorin		Terbutaline
	Furosemide	Methotrexate		Theophylline
	Methicillin	Mitoxantrone		
	Nafcillin	Vinblastine		
	Oxacillin			
	Pentamidine			
	Tazobactam			
	Clavulinate K			
	Tobramcine			
	Vancomycin			
	Imipenem			

Notices:

- Specific drugs or regimens referred to herein are strictly illustrative, and do not constitute claims by Ace Medical of effectiveness or recommended use.
 - Before use, consult the pharmaceutical manufacturer's guidelines for the medications or fluids being used with the delivery system.
 - As with all aseptically prepared products, the microbiological integrity of the product must be considered when assigning shelf life.
 - For complete instructions on proper use of the AutoFuser, refer to the device operating instructions supplied with the product.
-