

Spikes & Filters



RoweSpike I eco SV

Art.-No.	Material (body)	Filling volume (ml)	Connection	Tip	Hydrophilic filter	Specifics	PU ((pcs./box) / pal.)
A-6413	ABS - w	0.25	fLL	Spike	-	Safety Valve, green	300 / 16,800
A-6414	ABS - w	0.25	fLL	Spike	5 µm	Safety Valve, blue	300 / 16,800
A-6415	MABS - t	0.25	fLL	Spike	5 µm, chemo	Safety Valve, red	150 / 16,800
A-6416	MABS - t	0.25	fLL	Spike	chemo, free flow	Safety Valve, green	150 / 16,800

Made in Germany

ISO 13485



- For repeated withdrawal of solutions, drugs and toxic substances and ventilation of vial and containers
- Simple one-hand handling of safety snap cap
- Protection against contamination through safety snap cap
- Easy handling due to ergonomic design
- Low priming volume
- Low piercing forces
- Ventilation filter 0.1 µm PTFE membrane (hydrophobic bacteria filter)

w = white
t = transparent

Deep-drawn fluid opening to empty the bottle completely

